



DIGITAL BABY WEIGHING SCALE

- Micro-computer based system for ease in operation & accuracy
- Kg. to Pound conversion facility
- High bright red LED display for strain free reading
- Feather touch keys
- Hold with Acti-weigh technology to eliminate fluctuation in reading caused due to active baby
- Zero setting facility
- User friendly calibration facility
- Operating Voltage : 230 V AC +/- 10%, 50 Hz

Optional
Battery Backup


Model : 20	Model : 05
Capacity : 20 Kg.	Capacity : 5 Kg.
Resolution : 5 gms.	Resolution : 1 gm.
Accuracy : +/- 5 gms.	Accuracy : +/- 1 gm.

DIGITAL ADULT WEIGHING SCALE

- Micro-computer based system for ease in operation & accuracy
- Kg. to Pound conversion facility
- High bright red LED display for strain free reading
- Feather touch keys
- Hold with Acti-weigh technology to eliminate fluctuation in reading caused due to activeness
- Zero setting facility
- User friendly calibration facility
- Sturdy platform to weigh Adults
- Separate LED display module to be kept on table or wall mount
- Operating Voltage : 230 V AC +/- 10%, 50 Hz

Optional
Battery Backup

Model : 200
Capacity : 200 Kg.
Resolution : 20 gms.
Accuracy : +/- 20 gms.



Advance Mediplus Pvt. Ltd.

Showroom : M-72, Street No. 5, Shastri Nagar, New Delhi-110 052

Ph.: +91-11-23646624, 23647624

E-mail : contact@advancemedi.com Website : www.advancemedi.com

An ISO 9001:2008 certified company

