

3-IN-1

DIGITAL WEIGHING SCALE

DIGITAL WEIGHING SCALE

introduces Digital weighing scales which ensures accuracy in weight during active conditions caused due to movements.

- Micro-computer based system for ease in operation & accuracy
- Kg. to Pound conversion facility
- Hold with Acti-weight technology to eliminate fluctuation in reading caused due to active baby
- High bright red LED display for strain free reading
- Folding seat to weigh sitting baby
- Feather touch keys
- Zero setting facility
- Can weigh for Neonate, Infant & Adults
- Sturdy platform to weigh Adults
- Separate LED display module to be kept on table or wall mount
- User friendly calibration facility
- Operating Voltage : 230 V AC +/- 10%, 50 Hz

Optional
Battery Backup

Model : **3-in-1**
Capacity : 200 Kg.
Resolution : 20 gms.
Accuracy : +/- 20gms.



A unique system which fulfills all requirements of human weight by weighing in supine, sitting & standing position.

