

SIEMENS 300 & 300A RESPIRATORY VENTILATOR

PRODUCT DESCRIPTION

Refurbished Siemens 300 & 300A Servo Ventilators offer a wide range of ventilation modes together with an Automode function for greater treatment possibilities, shorter weaning times, and increased patient comfort.

Features:

- Neonatal to adult patient range
- Wide range of ventilation modes
- Controlled, Support, or Combined ventilation modes
- Easy to use
- Optional graphical display, and compressor
- Refurbished

Configurable Model:

- Siemens 300A Servo Ventilator
- Automode for fully adaptive patient interaction
- Neonate to adult patient range
- Wide range of ventilation modes: Controlled, Supported, Combined
- Easy to use
- Refurbished

Optional Accessories for Product:

- Carts
- Breathing Systems
- Compressors
- External Power Supply
- Humidifiers

SPECIFICATIONS

Dimensions

- Patient unit: 9.5' (W) x 14.6" (D) x 9.4' (H)
- User interface: 17.0" (W) x 15.9" (D) x 12.8" (H)

Weight: Approx. 52.9 lbs

Ventilation Modes:

- Controlled
- Supported
- Combined

Patient Range:

- Large/Medium/Small Animal

Flow Range: 0.1 ml/s - 3 l/s

Power Supply:

100, 120, 220, and 240 VAC \pm 10%, 50-60Hz

Alarms:

- Airway Pressure (Upper): 15-120 cm H₂O
- High Continuous Pressure: Set PEEP-level + 15 cm H₂O for more than 15 seconds.
- O₂ concentration: Set value \pm 6%. Lower alarm limit cannot go below 18% or exceed 90%.
- Expired Minute Volume (Upper Alarm Limit): Adult/Pediatric Range: 0-60 l/min, Neonate Range: 0-6 l/min.
- Expired Minute Volume (Lower Alarm Limit): Adult/Pediatric Range: 0.3-40 l/min, Neonate Range: 0.06-4 l/min.
- Apnea: Adult: 20 s. Pediatric: 15 s. Neonates: 10 s.

