

Infant T-Piece Resuscitator



The Infant T-Piece Resuscitator is a manually operated, gas-powered resuscitator designed to provide breaths at a set flow with consistent Peak Inspiratory Pressure (PIP) and Positive End Expiratory Pressure (PEEP) to infants.

The Infant T-Piece Resuscitator is an easy-to-use, gas-powered resuscitator, designed to provide safe, consistent and optimal resuscitation for infants. Recently it has been updated to further enhance functionality and usability while still allowing for breaths to be delivered at a set flow with consistent Peak Inspiratory Pressure (PIP) and Positive End Expiratory Pressure (PEEP).